DILATOMETRY 6 – 15



DILATOMETER CUTTING TOOL SPECIMEN PRECISION DCT

DESCRIPTION

This instrument, especially developed, cut by means of a special device at fixed measure of 25 mm or 50 mm the specimen after a first rough out to insert into the dilatometer. (The angles of the sample opposite faces must be properly rounded before

the use).Other sample measures on demand.Different from the standard

EQUIPMENT

- Lower tank and special devices made stainless steel AISI-304
- · Motor support in aluminium epoxi painted
- Diamond disc Ø80 mm
- Fixed speed 1000 rpm
- Vertical stroke 50 mm
- Cooling system by means upper tank and ball valve
- Safety micro-switch
- · Plexiglas protection
- Supply: Accumulator NI-MH 11,4 Vdc

ACCESSORIES AND SPARE PARTS

GT0953 Diamond disc Ø 80 mm